

For Research Use Only

Mouse 2B4/CD244/SLAMF4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84628-1-PBS



Basic Information

Catalog Number:

84628-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_018729.2

GeneID (NCBI):

18106

UNIPROT ID:

Q07763-1

Full Name:

CD244 natural killer cell receptor 2B4

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

242044C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01438-1, MOUSE CD244 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84628-2-PBS. Detection antibody: 84628-1-PBS. Standard: Eg1583. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP01438-2, MOUSE CD244 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84628-3-PBS. Detection antibody: 84628-1-PBS. Standard: Eg1583. Range: 0.781-50 ng/mL