

For Research Use Only

Mouse 2B4/CD244/SLAMF4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84628-3-PBS



Basic Information

Catalog Number: 84628-3-PBS	GenBank Accession Number: NM_018729.2	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 18106	CloneNo.: 242044G12
Source: Rabbit	UNIPROT ID: Q07763-1	
Isotype: IgG	Full Name: CD244 natural killer cell receptor 2B4	
Immunogen Catalog Number: EG1583	Calculated MW: 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01438-2, MOUSE CD244 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84628-3-PBS. Detection antibody: 84628-1-PBS. Standard: Eg1583. Range: 0.781-50 ng/mL.