

For Research Use Only

C12orf51 Recombinant antibody, PBS Only

Catalog Number: 84641-1-PBS



Basic Information

Catalog Number: 84641-1-PBS	GenBank Accession Number: BC006270	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 283450	CloneNo.: 241815A12
Source: Rabbit	Full Name: chromosome 12 open reading frame 51	
Isotype: IgG	Calculated MW: 3996 aa, 439 kDa	
Immunogen Catalog Number: AG14614	Observed MW: 500 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

The HECT domain E3 ubiquitin protein ligase 4 (HECTD4), is also known as C12orf51. It may be involved in protein modification or ubiquitination (PMID: 26891264). In addition, genetic polymorphisms of C12orf51 are related to alcohol consumption and type 2 diabetes (T2D) (PMID: 38523829).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84641-1-RR (C12orf51 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84641-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84641-1-RR against Human C12orf51 were performed. The affinity constant is below 1 pM.