

For Research Use Only

EPCR/CD201 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84648-2-PBS



Basic Information

Catalog Number: 84648-2-PBS	GenBank Accession Number: BC014451	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10544	CloneNo.: 241806D1
Source: Rabbit	UNIPROT ID: Q9UNN8	
Isotype: IgG	Full Name: protein C receptor, endothelial (EPCR)	
	Calculated MW: 238 aa, 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01480-2, PROCR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84648-2-PBS. Detection antibody: 84648-4-PBS. Standard: Eg2339. Range: 0.625-80 ng/mL