

For Research Use Only

Ubiquilin 1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84649-3-PBS



Basic Information

Catalog Number:

84649-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17530

GenBank Accession Number:

BC010066

GeneID (NCBI):

29979

UNIPROT ID:

Q9UMX0

Full Name:

ubiquilin 1

Calculated MW:

589 aa, 63 kDa

Purification Method:

Protein A purification

CloneNo.:

242092F12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01462-2, Ubiquitin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84649-2-PBS. Detection antibody: 84649-3-PBS. Standard: Ag17530. Range: 0.469-30 ng/mL.