

For Research Use Only

NEK4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84654-3-PBS



Basic Information

Catalog Number:

84654-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34929

GenBank Accession Number:

BC063044

GeneID (NCBI):

6787

UNIPROT ID:

P51957

Full Name:

NIMA (never in mitosis gene a)-related kinase 4

Calculated MW:

781 aa, 88 kDa

Purification Method:

Protein A purification

CloneNo.:

242081F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01486-1, NEK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84654-2-PBS. Detection antibody: 84654-3-PBS. Standard: Ag34929. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01486-2, NEK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84654-1-PBS. Detection antibody: 84654-3-PBS. Standard: Ag34929. Range: 0.625-80 ng/mL