

For Research Use Only

# PSPC1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84657-1-PBS



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
84657-1-PBS	BC014184	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	55269	242158B8
Source:	UNIPROT ID:	
Rabbit	Q8WXF1	
Isotype:	Full Name:	
IgG	paraspeckle component 1	
Immunogen Catalog Number:	Calculated MW:	
AG10136	523 aa, 59 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01461-1, PSPC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84657-1-PBS. Detection antibody: 84657-2-PBS. Standard:Ag10136. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01461-2, PSPC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84657-1-PBS. Detection antibody: 84657-3-PBS. Standard:Ag10136. Range: 0.781-100 ng/mL.