

For Research Use Only

# TRIM26 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84661-3-PBS



## Basic Information

Catalog Number:

84661-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25546

GenBank Accession Number:

BC024039

GeneID (NCBI):

7726

UNIPROT ID:

Q12899

Full Name:

tripartite motif-containing 26

Calculated MW:

62 kDa

Purification Method:

Protein A purification

CloneNo.:

242154C11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01471-2, TRIM26 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84661-1-PBS. Detection antibody: 84661-3-PBS. Standard: Ag25546. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01471-1, TRIM26 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84661-2-PBS. Detection antibody: 84661-3-PBS. Standard: Ag25546. Range: 0.313-40 ng/mL