

For Research Use Only

# COMMD4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84664-2-PBS



## Basic Information

Catalog Number:

84664-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34904

GenBank Accession Number:

BC068998

GeneID (NCBI):

54939

UNIPROT ID:

Q9HOA8

Full Name:

COMM domain containing 4

Purification Method:

Protein A purification

CloneNo.:

242175E3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01475-1, COMMD4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84664-3-PBS. Detection antibody: 84664-2-PBS. Standard: Ag34904. Range: 0.781-100 ng/mL