

For Research Use Only

C1orf123 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84665-2-PBS



Basic Information

Catalog Number:

84665-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36156

GenBank Accession Number:

BC010908

GeneID (NCBI):

54987

UNIPROT ID:

Q9NWX4

Full Name:

chromosome 1 open reading frame

123

Purification Method:

Protein A purification

CloneNo.:

242165B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean C1orf123 concentration was determined to be 579.7 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP01473-1, C1orf123 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84665-3-PBS. Detection antibody: 84665-2-PBS. Standard: Ag36156. Range: 0.781-100 ng/mL