

For Research Use Only

# FGF13 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84666-3-PBS



## Basic Information

Catalog Number: 84666-3-PBS	GenBank Accession Number: BC034340	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2258	CloneNo.: 242153D9
Source: Rabbit	UNIPROT ID: Q92913	
Isotype: IgG	Full Name: fibroblast growth factor 13	
Immunogen Catalog Number: AG23965	Calculated MW: 28 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01481-2, FGF13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84666-1-PBS. Detection antibody: 84666-3-PBS. Standard: Ag23965. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01481-1, FGF13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84666-2-PBS. Detection antibody: 84666-3-PBS. Standard: Ag23965. Range: 0.781-100 ng/mL.