

For Research Use Only

RPL11 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84670-1-PBS



Basic Information

Catalog Number:

84670-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9356

GenBank Accession Number:

BC018970

GeneID (NCBI):

6135

UNIPROT ID:

P62913

Full Name:

ribosomal protein L11

Calculated MW:

177 aa, 20 kDa

Purification Method:

Protein A purification

CloneNo.:

242147F7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01493-1, RPL11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84670-1-PBS. Detection antibody: 84670-2-PBS. Standard: Ag9356. Range: 0.625-80 ng/mL