

For Research Use Only

DIAPH2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84671-2-PBS



Basic Information

Catalog Number:

84671-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35280

GenBank Accession Number:

NM_006729.4

GeneID (NCBI):

1730

UNIPROT ID:

O60879

Full Name:

diaphanous homolog 2 (Drosophila)

Calculated MW:

126kDa

Purification Method:

Protein A purification

CloneNo.:

242176B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01492-2, DIAPH2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84671-1-PBS. Detection antibody: 84671-2-PBS. Standard: Ag35280. Range: 0.156-20 ng/mL