

For Research Use Only

# CST5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84678-1-PBS



## Basic Information

Catalog Number:

84678-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34412

GenBank Accession Number:

BC069514

GeneID (NCBI):

1473

UNIPROT ID:

P28325

Full Name:

cystatin D

Purification Method:

Protein A purification

CloneNo.:

241564C2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01490-2, CST5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84678-1-PBS. Detection antibody: 84678-2-PBS. Standard: Ag34412. Range: 0.781-100 ng/mL.