

For Research Use Only

CDC40 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84679-2-PBS



Basic Information

Catalog Number:

84679-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36151

GenBank Accession Number:

BC126114

GeneID (NCBI):

51362

UNIPROT ID:

O60508

Full Name:

cell division cycle 40 homolog (S. cerevisiae)

Purification Method:

Protein A purification

CloneNo.:

242034H4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01501-1, CDC40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84679-2-PBS. Detection antibody: 84679-1-PBS. Standard: Ag36151. Range: 0.313-40 ng/mL