

For Research Use Only

# FOXM1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84681-3-PBS



## Basic Information

Catalog Number:

84681-3-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3729

GenBank Accession Number:

BC006192

GeneID (NCBI):

2305

UNIPROT ID:

Q08050

Full Name:

forkhead box M 1

Calculated MW:

83 kDa

Purification Method:

Protein A purification

CloneNo.:

242077D9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01500-2, FOXM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84681-1-PBS. Detection antibody: 84681-3-PBS. Standard: Ag3729. Range: 0.625-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 84681-3-PBS against Human FOXM1 was performed. The affinity constant is 0.712 nM.