

For Research Use Only

Semaphorin 7A Recombinant antibody, PBS Only (Capture)

Catalog Number: 84682-2-PBS



Basic Information

Catalog Number:

84682-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003612.5

GeneID (NCBI):

8482

UNIPROT ID:

O75326-1

Full Name:

semaphorin 7A, GPI membrane
anchor (John Milton Hagen blood
group)

Calculated MW:

75kDa

Purification Method:

Protein A purification

CloneNo.:

241814B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01498-1, SEMA7A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84682-2-PBS. Detection antibody: 84682-1-PBS. Standard: Eg2005. Range: 0.781-100 ng/mL.