

For Research Use Only

Siglec-9 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84683-3-PBS



Basic Information

Catalog Number:

84683-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014441.3

GeneID (NCBI):

27180

UNIPROT ID:

Q9Y336-1

Full Name:

sialic acid binding Ig-like lectin 9

Calculated MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

241799G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01502-1, Siglec-9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84683-3-PBS. Detection antibody: 84683-2-PBS. Standard: Eg1084. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01502-2, Siglec-9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84683-3-PBS. Detection antibody: 84683-1-PBS. Standard: Eg1084. Range: 0.156-20 ng/mL