

For Research Use Only

# CCDC111 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84704-3-PBS



## Basic Information

Catalog Number:

84704-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31478

GenBank Accession Number:

BC064600

GeneID (NCBI):

201973

UNIPROT ID:

Q96LW4

Full Name:

coiled-coil domain containing 111

Purification Method:

Protein A purification

CloneNo.:

242119G11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01510-1, CCDC111 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84704-3-PBS. Detection antibody: 84704-1-PBS. Standard: Ag31478. Range: 0.078-10 ng/mL.



Cytometric bead array standard curve of MP01510-2, CCDC111 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84704-3-PBS. Detection antibody: 84704-2-PBS. Standard: Ag31478. Range: 0.078-10 ng/mL.