

For Research Use Only

TAF1D Recombinant antibody, PBS Only (Detector)

Catalog Number: 84712-1-PBS



Basic Information

Catalog Number: 84712-1-PBS	GenBank Accession Number: BC001972	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 79101	CloneNo.: 242087C9
Source: Rabbit	UNIPROT ID: Q9H5J8	
Isotype: IgG	Full Name: TATA box binding protein (TBP)- associated factor, RNA polymerase I, D, 41kDa	
Immunogen Catalog Number: AG10014	Calculated MW: 32 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01514-1, TAF1D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84712-3-PBS. Detection antibody: 84712-1-PBS. Standard: Ag10014. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01514-2, TAF1D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84712-2-PBS. Detection antibody: 84712-1-PBS. Standard: Ag10014. Range: 0.625-80 ng/mL