

For Research Use Only

WASP Recombinant antibody, PBS Only (Capture)

Catalog Number: 84714-1-PBS



Basic Information

Catalog Number:

84714-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1444

GenBank Accession Number:

BC012738

GeneID (NCBI):

7454

UNIPROT ID:

P42768

Full Name:

Wiskott-Aldrich syndrome (eczema-thrombocytopenia)

Calculated MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

242074A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01520-1, WASP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84714-1-PBS. Detection antibody: 84714-2-PBS. Standard: Ag1444. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01520-2, WASP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84714-1-PBS. Detection antibody: 84714-3-PBS. Standard: Ag1444. Range: 0.625-80 ng/mL