

For Research Use Only

# Rat Ferritin light chain 1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84731-2-PBS



## Basic Information

Catalog Number:

84731-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_022500.4

GeneID (NCBI):

29292

UNIPROT ID:

P02793

Full Name:

ferritin, light polypeptide

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

242299B7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01511-1, RAT Ferritin light chain 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84731-1-PBS. Detection antibody: 84731-2-PBS. Standard: Eg2434. Range: 0.625-80 ng/mL