

For Research Use Only

Rat Ferritin light chain 1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84731-3-PBS



Basic Information

Catalog Number:

84731-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_022500.4

GeneID (NCBI):

29292

UNIPROT ID:

P02793

Full Name:

ferritin, light polypeptide

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

242299E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01511-2, RAT Ferritin light chain 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84731-3-PBS. Detection antibody: 84731-4-PBS. Standard: Eg2434. Range: 0.625-80 ng/mL