

For Research Use Only

# EPB41L4B Recombinant antibody, PBS Only (Detector)

Catalog Number: 84732-1-PBS



## Basic Information

Catalog Number:

84732-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36353

GenBank Accession Number:

NM\_00138562

GeneID (NCBI):

54566

UNIPROT ID:

Q9H329

Full Name:

erythrocyte membrane protein band

4.1 like 4B

Calculated MW:

99kDa

Purification Method:

Protein A purification

CloneNo.:

242186D2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01525-1, EPB41L4B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84732-2-PBS. Detection antibody: 84732-1-PBS. Standard: Ag36353. Range: 0.781-100 ng/mL