

For Research Use Only

IDI2 Recombinant antibody, PBS Only

Catalog Number: 84746-1-PBS



Basic Information

Catalog Number:

84746-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10098

GenBank Accession Number:

BC017778

GeneID (NCBI):

91734

UNIPROT ID:

Q9BXS1

Full Name:

isopentenyl-diphosphate delta
isomerase 2

Calculated MW:

227 aa, 27 kDa

Observed MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

241885H3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84746-1-RR (IDI2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84746-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84746-1-RR against Human IDI2 were performed. The affinity constant is 6.67 nM.