

For Research Use Only

# KLRB1/CD161 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84755-3-PBS



## Basic Information

Catalog Number:

84755-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002258.3

GeneID (NCBI):

3820

ENSEMBL Gene ID:

ENSG00000111796

UNIPROT ID:

Q12918

Full Name:

killer cell lectin-like receptor  
subfamily B, member 1

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

242140F10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01539-1, KLRB1/CD161 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84755-1-PBS. Detection antibody: 84755-3-PBS. Standard: Eg1788. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01539-2, KLRB1/CD161 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84755-2-PBS. Detection antibody: 84755-3-PBS. Standard: Eg1788. Range: 0.156-20 ng/mL