

For Research Use Only

TXNDC5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84756-2-PBS



Basic Information

Catalog Number:

84756-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13965

GenBank Accession Number:

BC001199

GeneID (NCBI):

81567

UNIPROT ID:

Q8NBS9

Full Name:

thioredoxin domain containing 5
(endoplasmic reticulum)

Calculated MW:

432 aa, 48 kDa

Purification Method:

Protein A purification

CloneNo.:

242195A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01536-2, TXNDC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84756-3-PBS. Detection antibody: 84756-2-PBS. Standard: Ag13965. Range: 1.56-100 ng/mL.