

For Research Use Only

CAB39 Recombinant antibody, PBS Only

Catalog Number: 84760-4-PBS



Basic Information

Catalog Number:

84760-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36153

GenBank Accession Number:

BC020570

GeneID (NCBI):

51719

UNIPROT ID:

Q9Y376

Full Name:

calcium binding protein 39

Observed MW:

39-40 kDa

Purification Method:

Protein A purification

CloneNo.:

242094G11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

CAB39 (Calcium-binding protein 39), also called MO25, is a component of the complex that binds and activates STK11/LKB1. It is required to stabilize the interaction between CAB39/MO25 (CAB39/MO25alpha or CAB39L/MO25beta) and STK11/LKB1 (PMID: 19892943).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84760-4-RR (CAB39 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84760-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84760-4-RR against Human CAB39 were performed. The affinity constant is 0.304 nM.