

For Research Use Only

SPTAN1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84761-1-PBS



Basic Information

Catalog Number:

84761-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36472

GenBank Accession Number:

NM_003127.3

GeneID (NCBI):

6709

UNIPROT ID:

Q13813

Full Name:

spectrin, alpha, non-erythrocytic 1
(alpha-fodrin)

Calculated MW:

285kDa, 2472aa

Purification Method:

Protein A purification

CloneNo.:

242104B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01560-1, SPTAN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84761-2-PBS. Detection antibody: 84761-1-PBS. Standard: Ag36472. Range: 0.469-60 ng/mL.