

For Research Use Only

PDIA2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84765-3-PBS



Basic Information

Catalog Number:

84765-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35301

GenBank Accession Number:

NM_006849.2

GeneID (NCBI):

64714

UNIPROT ID:

Q13087

Full Name:

protein disulfide isomerase family A, member 2

Calculated MW:

58 kDa

Purification Method:

Protein A purification

CloneNo.:

242294H6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01546-1, PDIA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84765-3-PBS. Detection antibody: 84765-1-PBS. Standard: Ag35301. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01546-2, PDIA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84765-3-PBS. Detection antibody: 84765-2-PBS. Standard: Ag35301. Range: 0.156-20 ng/mL