

For Research Use Only

# HMGCL Recombinant antibody, PBS Only (Capture)

Catalog Number: 84769-2-PBS



## Basic Information

Catalog Number:

84769-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10426

GenBank Accession Number:

BC010570

GeneID (NCBI):

3155

UNIPROT ID:

P35914

Full Name:

3-hydroxymethyl-3-methylglutaryl-Coenzyme A lyase

Calculated MW:

325 aa, 34 kDa

Purification Method:

Protein A purification

CloneNo.:

242224H1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01555-2, HMGCL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84769-2-PBS. Detection antibody: 84769-1-PBS. Standard: Ag10426. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01555-1, HMGCL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84769-2-PBS. Detection antibody: 84769-3-PBS. Standard: Ag10426. Range: 0.313-40 ng/mL