

For Research Use Only

# HNMT Recombinant antibody, PBS Only (Capture)

Catalog Number: 84777-1-PBS



## Basic Information

Catalog Number:

84777-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2462

GenBank Accession Number:

BC020677

GeneID (NCBI):

3176

UNIPROT ID:

P50135

Full Name:

histamine N-methyltransferase

Calculated MW:

292 aa, 33 kDa

Purification Method:

Protein A purification

CloneNo.:

242159D11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01567-1, HNMT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84777-1-PBS. Detection antibody: 84777-2-PBS. Standard: Ag2462. Range: 0.781-100 ng/mL.