

For Research Use Only

# clAP2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84779-3-PBS



## Basic Information

Catalog Number:

84779-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20500

GenBank Accession Number:

BC027485

GeneID (NCBI):

330

UNIPROT ID:

Q13489

Full Name:

baculoviral IAP repeat-containing 3

Calculated MW:

604 aa, 68 kDa

Purification Method:

Protein A purification

CloneNo.:

242214G1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01556-1, BIRC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84779-3-PBS. Detection antibody: 84779-1-PBS. Standard: Ag20500. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP01556-2, BIRC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84779-3-PBS. Detection antibody: 84779-2-PBS. Standard: Ag20500. Range: 0.313-20 ng/mL