

For Research Use Only

MDM2/HDM2 Recombinant antibody, PBS Only

Catalog Number: 84780-3-PBS



Basic Information

Catalog Number:

84780-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27255

GenBank Accession Number:

NM_001145337.2

GeneID (NCBI):

4193

UNIPROT ID:

Q00987

Full Name:

Mdm2 p53 binding protein homolog
(mouse)

Calculated MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

242080D3

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

