

For Research Use Only

# CNTLN Recombinant antibody, PBS Only (Detector)

Catalog Number: 84781-3-PBS



## Basic Information

Catalog Number:

84781-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36354

GenBank Accession Number:

GeneID (NCBI):

54875

UNIPROT ID:

Q9NXG0

Full Name:

centlein, centrosomal protein

Calculated MW:

162kDa

Purification Method:

Protein A purification

CloneNo.:

242187A4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01557-2, CNTLN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84781-2-PBS. Detection antibody: 84781-3-PBS. Standard: Ag36354. Range: 0.313-40 ng/mL