

For Research Use Only

Mouse APRIL/TNFSF13 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 84789-1-PBS



Basic Information

Catalog Number: 84789-1-PBS	GenBank Accession Number: NM_023517.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 69583	CloneNo.: 241666A4
Isotype: IgG	UNIPROT ID: Q9D777	
Immunogen Catalog Number: EG1561	Full Name: tumor necrosis factor (ligand) superfamily, member 13	
	Calculated MW: 27kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01586-1, MOUSE APRIL/TNFSF13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84789-1-PBS. Detection antibody: 84789-2-PBS. Standard: Eg1561. Range: 0.781-25 ng/mL