

For Research Use Only

Mouse Siglec-12 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84790-1-PBS



Basic Information

Catalog Number:

84790-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_031181.2

GeneID (NCBI):

83382

UNIPROT ID:

Q91Y57

Full Name:

sialic acid binding Ig-like lectin E

Calculated MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

242248A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01535-1, MOUSE Siglec-12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84790-2-PBS. Detection antibody: 84790-1-PBS. Standard: Eg2733. Range: 3.125-100 ng/mL.