

For Research Use Only

Neurophysin 1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 84791-3-PBS



Basic Information

Catalog Number:

84791-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2237

GenBank Accession Number:

BC101841

GeneID (NCBI):

5020

UNIPROT ID:

P01178

Full Name:

oxytocin, prepropeptide

Calculated MW:

125 aa, 13 kDa

Purification Method:

Protein A purification

CloneNo.:

242251F8

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

