

For Research Use Only

Mouse gp130/IL6ST Recombinant antibody, PBS Only (Capture)

Catalog Number: 84794-3-PBS



Basic Information

Catalog Number: 84794-3-PBS	GenBank Accession Number: NM_010560.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 16195	CloneNo.: 242236D9
Source: Rabbit	UNIPROT ID: Q00560	
Isotype: IgG	Full Name: interleukin 6 signal transducer	
	Calculated MW: 102kDa, 917aa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01592-1, MOUSE gp130/IL6ST Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84794-3-PBS. Detection antibody: 84794-1-PBS. Standard: Eg1631. Range: 0.469-30 ng/mL