

For Research Use Only

JMJD5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84804-3-PBS



Basic Information

Catalog Number:

84804-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15115

GenBank Accession Number:

BC027911

GeneID (NCBI):

79831

UNIPROT ID:

Q8N371

Full Name:

jumonji domain containing 5

Calculated MW:

416 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

242097H4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01571-2, JMJD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84804-3-PBS. Detection antibody: 84804-2-PBS. Standard: Ag15115. Range: 0.625-80 ng/mL