

For Research Use Only

SIRT7 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84805-2-PBS



Basic Information

Catalog Number:

84805-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31794

GenBank Accession Number:

NM_016538

GeneID (NCBI):

51547

UNIPROT ID:

Q9NRC8

Full Name:

sirtuin (silent mating type
information regulation 2 homolog) 7
(S. cerevisiae)

Calculated MW:

45KD

Purification Method:

Protein A purification

CloneNo.:

242167A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01572-1, SIRT7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84805-3-PBS. Detection antibody: 84805-2-PBS. Standard: Ag31794. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP01572-2, SIRT7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84805-1-PBS. Detection antibody: 84805-2-PBS. Standard: Ag31794. Range: 3.125-100 ng/mL.