

For Research Use Only

Mouse Fc-epsilon RI-alpha (FcERI) Recombinant antibody, PBS Only

Catalog Number: 84816-1-PBS



Basic Information

Catalog Number:

84816-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1740

GenBank Accession Number:

NM_010184.2

GeneID (NCBI):

14125

UNIPROT ID:

P20489

Full Name:

Fc receptor, IgE, high affinity I, alpha polypeptide

Calculated MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

242008A3

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

