

For Research Use Only

Mouse CD229/SLAMF3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84817-1-PBS



Basic Information

Catalog Number:

84817-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1586

GenBank Accession Number:

NM_001277968.1

GeneID (NCBI):

17085

UNIPROT ID:

Q01965

Full Name:

lymphocyte antigen 9

Calculated MW:

73 kDa

Purification Method:

Protein A purification

CloneNo.:

242312A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01599-2, MOUSE CD229/SLAMF3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84817-1-PBS. Detection antibody: 84817-2-PBS. Standard: Eg1586. Range: 1..563-100 ng/mL