

For Research Use Only

Mouse CD229/SLAMF3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84817-2-PBS



Basic Information

Catalog Number:

84817-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001277968.1

GeneID (NCBI):

17085

UNIPROT ID:

Q01965

Full Name:

lymphocyte antigen 9

Calculated MW:

73kDa

Purification Method:

Protein A purification

CloneNo.:

242312D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01599-1, MOUSE CD229/SLAMF3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84817-3-PBS. Detection antibody: 84817-2-PBS. Standard: Eg1586. Range: 0.781-25 ng/mL



Cytometric bead array standard curve of MP01599-2, MOUSE CD229/SLAMF3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84817-1-PBS. Detection antibody: 84817-2-PBS. Standard: Eg1586. Range: 1..563-100 ng/mL