

For Research Use Only

# Aggrecan Recombinant antibody, PBS Only (Detector)

Catalog Number: 84819-2-PBS



## Basic Information

Catalog Number: 84819-2-PBS	GenBank Accession Number: BC036445	Purification Method: Protein A purification
Concentration: 1000 µg/ml	GeneID (NCBI): 176	CloneNo.: 242281E10
Source: Rabbit	UNIPROT ID: P16112	
Isotype: IgG	Full Name: aggrecan	
Immunogen Catalog Number: AG4985	Calculated MW: 250 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01577-1, Aggrecan Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84819-3-PBS. Detection antibody: 84819-2-PBS. Standard: Ag4985. Range: 0.313-40 ng/mL