

For Research Use Only

Mouse Ly6g Recombinant antibody, PBS Only (Capture)

Catalog Number:84820-1-PBS



Basic Information

Catalog Number:

84820-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001310438.1

GeneID (NCBI):

546644

UNIPROT ID:

P35461

Full Name:

lymphocyte antigen 6 complex, locus G

Calculated MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

242141B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01602-2, MOUSE Ly-6G Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84820-1-PBS. Detection antibody: 84820-2-PBS. Standard: Eg2609. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP01602-1, MOUSE Ly-6G Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84820-1-PBS. Detection antibody: 84820-3-PBS. Standard: Eg2609. Range: 0.625-80 ng/mL.