

For Research Use Only

# EZH2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84824-5-PBS



## Basic Information

<b>Catalog Number:</b> 84824-5-PBS	<b>GenBank Accession Number:</b> BC010858	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 2146	<b>CloneNo.:</b> 242309E8
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q15910	
<b>Isotype:</b> IgG	<b>Full Name:</b> enhancer of zeste homolog 2 (Drosophila)	
<b>Immunogen Catalog Number:</b> AG16789	<b>Calculated MW:</b> 751 aa, 86 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP01605-3, Human EZH2 Recombinant Matched Antibody Pair - PBS only. 84824-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag16789. 84824-5-PBS was HRP conjugated as the detection antibody. Range: 0.938-60 ng/mL

The mean EZH2 concentration was determined to be 412.79 ng/mL in Raji cell extract based on a 1.0 mg/mL extract load.