

For Research Use Only

# DGKZ Recombinant antibody, PBS Only (Capture)

Catalog Number: 84830-2-PBS



## Basic Information

Catalog Number:

84830-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37336

GenBank Accession Number:

NM\_001199267.1

GeneID (NCBI):

8525

UNIPROT ID:

Q13574

Full Name:

diacylglycerol kinase, zeta 104kDa

Calculated MW:

104kDa, 928aa

Purification Method:

Protein A purification

CloneNo.:

241444C8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01606-2, DGKZ Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84830-2-PBS. Detection antibody: 84830-1-PBS. Standard: Ag37336. Range: 0.313-40 ng/mL