

For Research Use Only

DGKZ Recombinant antibody, PBS Only (Detector)

Catalog Number: 84830-4-PBS



Basic Information

Catalog Number: 84830-4-PBS	GenBank Accession Number: NM_001199267.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8525	CloneNo.: 241444G6
Source: Rabbit	UNIPROT ID: Q13574	
Isotype: IgG	Full Name: diacylglycerol kinase, zeta 104kDa	
Immunogen Catalog Number: AG37336	Calculated MW: 104kDa, 928aa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01606-1, DGKZ Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84830-3-PBS. Detection antibody: 84830-4-PBS. Standard: Ag37336. Range: 0.313-40 ng/mL