

For Research Use Only

REST Recombinant antibody, PBS Only (Capture)

Catalog Number: 84846-1-PBS



Basic Information

Catalog Number: 84846-1-PBS	GenBank Accession Number: BC132859	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5978	CloneNo.: 241789A9
Source: Rabbit	UNIPROT ID: Q13127	
Isotype: IgG	Full Name: RE1-silencing transcription factor	
Immunogen Catalog Number: AG17707	Calculated MW: 1097 aa, 122 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01629-1, REST Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84846-1-PBS. Detection antibody: 84846-2-PBS. Standard: Ag17707. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01629-2, REST Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84846-1-PBS. Detection antibody: 84846-3-PBS. Standard: Ag17707. Range: 0.156-20 ng/mL